RACK & TANK CEANER

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RACK & TANK CEANER

Rack & Tank Cleaner is a powdered special cleaner for developing machines in X-ray, Micro, Prepress and Black/White.

Rack & Tank Cleaner is **not suitable** for developer tanks with **hydroquinone-free developers**, likewise this cleaner must expressly **not** be used in **colour developer tanks**.

Pack for 4 x 10 L REF 17082 4 x Part A and Part B

PRE-CLEANING

Empty the tank to be cleaned and dispose of the used chemicals.

Fill the tank with water, switch on the machine briefly to activate the recirculation and then drain the water. Close the tank drain.

PREPARATION OF THE CLEANER

Avoid any formation of dust.

The cleaner should preferably be prepared outside the laboratory rooms to avoid possible contamination of the rooms with powder dust.

The preparation of the cleaner solution should not be carried out in the processor tank, but in an external mixer or in a separate container.

Dissolve Part A and Part B successively in approx. 9.6 litres of water to get approx. 10 litres of cleaning solution. Water temperature: max. 25°C.

Stir constantly until the powder is dissolved. If the tank to be cleaned is smaller or larger than 10 litres, the cleaner solution can be max. 30 % more concentrated or max. 20 % more diluted. Likewise, this applies to tanks with a multiple size of 10 litres.

APPLICATION

Insert the rack, remove the filter insert. Close the tank drain. Fill the cleaner to the overflow. Switch on the machine for ca. 30 minutes - max. 1 hour, not longer.

Then switch off the machine, empty the tank and collect the used cleaning solution for proper disposal. Close the tank drain.

Fill the tank with water, switch on the machine, let it run for 2-3 minutes, then switch it off, drain the water and close the tank drain. Repeat rinsing for 3 times.

Ready again for operation.